

# PULSE STREAMER 快速指南

本文旨在帮助客户初步了解 PULSE STREAMER 的使用。详细的参数和进一步的深入使用，请仔细查阅产品手册或联系相关工作人员。阅读本文前，请先查看《PULSE STREAMER 安装培训准备清单》并完成准备工作。在准备工作中，如果遇到任何问题，欢迎联系我们。



## 仪器连接

设备连接上电源和网线之后，设置网络连接：

# Example:

# using default hostname


```
ps = PulseStreamer('pulsestreamer')
```

# using fallback IP

```
ps = PulseStreamer('169.254.8.2')
```

## 软件操作

打开 Pulse Streamer 应用软件，双击设备 ID，



打开设置界面，如下图所示，在 STEADY STATE OUTPUT——FINAL STATE 中选择 Digital 0，在 Built-in 中的 DIGITAL 0 中选择 Pulse，点击上方的运行，此时 STATUS 中的指示灯将会变为 Streaming Output。最后，将仪器 Digital Out 0 接到示波器上观察输出的脉冲信号。

**PS** Pulse Streamer

File Help

**STATUS**

☒ Streaming output

**TRIGGER SETTINGS**

Trigger mode Immediate

Repeat indefinitely ☒

Repetition count 1

**STEADY STATE OUTPUT**

**FINAL STATE**

Digital ☒ 0 ☐ 1 ☐ 2 ☐ 3  
☐ 4 ☐ 5 ☐ 6 ☐ 7

Analog 0 0.0000 V

Analog 1 0.0000 V

**Built-in Arbitrary**



**DIGITAL 0**

Period 1,000 ns

Duration 500 ns

Time offset 0 ns

**DIGITAL 1**

State Constant

**DIGITAL 2**

State Constant

**DIGITAL 3**

State Constant

**DIGITAL 4**

State Constant

**DIGITAL 5**

State Constant

## 简单信号生成

以下程序表示从 digital output 0 输出无限循环的单个脉冲信号。

### Python

```
# import API classes into the current namespace
from pulsestreamer import PulseStreamer, Sequence

# A pulse with 10,ts HIGH and 30,ts LOW levels
pattern = [(10000, 1), (30000, 0)]

# Connect to Pulse Streamer
ip = 'pulsestreamer'

ps = PulseStreamer(ip)

# Create a sequence object
sequence = ps.createSequence()

# Create sequence and assign pattern to digital channel 0 sequence.setDigital(0, pattern)

# Stream the sequence and repeat it indefinitely
n_runs = PulseStreamer.REPEAT_INFINITY

ps.stream(sequence, n_runs )
```

## MATLAB

```
% import API classes into the current namespace
import PulseStreamer.*

% A pulse with 10,ts HIGH and 30,ts LOW levels
pattern = {10000, 1; 30000, 0};

% Connect to Pulse Streamer
ip = 'pulsestreamer';

ps = PulseStreamer(ip);

% Create a sequence object
sequence = ps.createSequence();

% Assign pulse pattern to digital channel 0
sequence.setDigital(0, pattern);

% Stream the sequence and repeat it indefinitely
n_runs = PulseStreamer.REPEAT_INFINITY; % endless streaming ps.stream(sequence, n_runs)
```

